



Our doping kits are used extensively in sporting events worldwide, including: Athletics, Basketball, Boxing, Cycling, Soccer, Golf, Ice Hockey, Judo, Luge, Marathon, Rugby, Sailing, Ski jumping, Skiing, Snowboarding, Summer Olympics, Swimming, Tennis, Triathlon, Weightlifting, Winter Olympics, and more...

**versapak**  
Doping Control Limited

## Secure & Tamper Evident Doping Control Products

Stringently tested by independent security experts, Versapak's unique tamper-evident containers provide a secure and user-friendly means with which doping control officers can collect and transport samples to laboratories worldwide.



Our anti-doping products include urine sampling kits (complete and partial), urine collection vessels and blood sampling kits, as well as security holdalls and security seals which facilitate the secure transportation of the kits.

Versapak Doping Control Ltd is dedicated to the design, manufacture and supply of doping control products. These products are used exclusively in the worldwide fight against the use of drugs in sport.

The Versapak Centre, Centurion Way, Erith, Kent, DA18 4AF, UK  
Versapak Doping Control Ltd is part of The Versapak Group  
[www.theversapakgroup.com](http://www.theversapakgroup.com)

Phone: +44 (0)20 8333 5300 | Fax: +44 (0)20 8311 0423  
Email: [info@versapakdopingcontrol.com](mailto:info@versapakdopingcontrol.com)  
Website: [www.versapakdopingcontrol.com](http://www.versapakdopingcontrol.com)

# Versapak Doping Control Ltd

## Anti Doping Product Details

### Urine Sampling Kit

#### DTC2

This Urine Sampling Kit has been tried and tested both by independent security experts and by anti doping officers around the world, who have used this kit since its creation in 1991.



#### Key components of DTC2:

##### 2 x screw capped 100ml capacity glass bottles

- Colour coded screw caps help to **identify A and B samples**
- **Tamper evident bottle labels** identify A and B sample numbers (which are also embossed on the security canister)
- When filled and secured inside the security canister, the glass bottle is not visible, creating an additional obstacle to tampering

##### 2 x HDPE security canisters

- The separate security canister provides an **extra level of security and protection during transport**
- Manufactured from high density polyethalene (HDPE) means that if breached, this **canister cannot be covertly repaired**
- The canister lid has an interference fit with a double bore seal - once closed, the **canister cannot be opened without showing visible evidence**
- Colour coded canisters **help to identify A and B samples**

##### 2 x absorbent pads

- In the unlikely event of the bottle breaking during transport, the absorbent pad will soak up any spilled urine

##### 1 x tamper evident outer bag

- Prevents tampering prior to sample collection
- Kit components can be easily viewed and identified prior to opening

# Partial Sampling Urine Kit BPK2



Manufactured with the same security protocols as the DTC2, this screw capped 100ml glass bottle with black lid and black HDPE security canister is used for securely storing and transporting partial urine samples.

### Urine Collection Vessel UCV



Used for the safe, spill-proof collection of the sample, and its transfer into the secure kit, this collection vessel has a push-fit lid which is opened by a ring pull. A calibration scale indicates sample volume, and a unique code embossed on each vessel ensures accurate identification.

### Security Holdall CCBX20

These purpose-built courier holdalls hold 20 sealed canisters. Made from heavy duty weatherproof PVC, they are tamper evident when sealed with a Versapak button seal. Designed for long life, they are reusable over 2,000 times.

#### Versapak Button Seal



Phone: +44 (0)20 8333 5300 | Fax: +44 (0)20 8311 0423

Email: [info@versapakdopingcontrol.com](mailto:info@versapakdopingcontrol.com)

Website: [www.versapakdopingcontrol.com](http://www.versapakdopingcontrol.com)